

丙型肝炎病毒 NS5A 激活蛋白 9 抗体

产品货号： mlR16961

英文名称： NS5ATP9

中文名称： 丙型肝炎病毒 NS5A 激活蛋白 9 抗体

别名： KIAA0101 HCV NS5A transactivated protein 9; HCV NS5A-transactivated protein 9; Hepatitis C virus NS5A transactivated protein 9; Hepatitis C virus NS5A-transactivated protein 9; KIAA0101; L5; NS5 ATP9; NS5ATP9; OEATC 1; OEATC; OEATC-1; OEATC1; Overexpressed in anaplastic thyroid carcinoma 1; p15(PAF); p15PAF; PAF; PAF_HUMAN; PCNA associated factor; PCNA-associated factor.

研究领域： 细胞生物 细胞凋亡 细菌及病毒 新陈代谢

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Pig, Cow, Horse, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 12kDa

细胞定位： 细胞核

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from human NS5ATP9:21-80/111

亚型： IgG

纯化方法： affinity purified by Protein A

储存液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.



PubMed : PubMed