

5 号染色体开放阅读框 42 抗体

产品货号： mlR15204

英文名称： C5ORF42

中文名称： 5 号染色体开放阅读框 42 抗体

别名： Chromosome 5 open reading frame 42; FLJ13231; Hypothetical protein LOC65250; Transmembrane protein ENSP00000382582; Uncharacterized protein C5orf42; CE042_HUMAN.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Cow, Sheep,

产品应用： ELISA=1:500-1000 Flow-Cyt=1 μ g/Test

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 362kDa

细胞定位： 细胞膜

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from human C5ORF42:1-100/3197

亚型： IgG

纯化方法： affinity purified by Protein A

储存液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed： PubMed

产品图片

