

二磷酸腺苷核糖基化因子 6 相互作用蛋白 抗体

产品货号： mIR10561

英文名称： ARL6

中文名称： 二磷酸腺苷核糖基化因子 6 相互作用蛋白抗体

别名： ADP ribosylation factor like 6; ADP ribosylation factor like protein 6; ADP-ribosylation factor-like protein 6; Arl6; ARL6_HUMAN; Bardet Biedl syndrome 3 protein; Bardet-Biedl syndrome 3 protein; BBS3; MGC32934.

研究领域： 细胞生物 信号转导 细胞凋亡 G 蛋白信号

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Horse,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500
(石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 20kDa

细胞定位： 细胞浆 细胞膜

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from human ARL6:31-120/186

亚型： IgG

纯化方法： affinity purified by Protein A

储存液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : The protein encoded by this gene belongs to the ARF family of GTP-binding proteins. ARF proteins are important regulators of cellular traffic and are the founding members of an expanding family of homologous proteins and genomic sequences. They depart from other small GTP-binding proteins by a unique structural device that implements front-back communication from the N-terminus to the nucleotide-binding site. Studies of the mouse ortholog of this protein suggest an involvement in protein transport, membrane trafficking, or cell signaling during hematopoietic maturation. Alternative splicing occurs at this locus and two transcript variants encoding the same protein have been described. [provided by RefSeq, Jul 2008].

Function:

Involved in membrane protein trafficking at the base of the ciliary organelle. Mediates recruitment onto plasma membrane of the BBSome complex which would constitute a coat complex required for sorting of specific membrane proteins to the primary cilia. May regulate cilia assembly and disassembly and subsequent ciliary signaling events such as the Wnt signaling cascade. Isoform 2 may be required for proper retinal function and organization.

Subunit:

Interacts with SEC61B, ARL6IP1, ARL6IP2, ARL6IP3, ARL6IP4 ARL6IP5 and ARL6IP6. Interacts (GTP-bound form) with the BBSome a complex that contains BBS1, BBS2, BBS4, BBS5, BBS7, BBS8, BBS9 and BBIP10.

Subcellular Location:

Cell projection, cilium membrane; Peripheral membrane protein; Cytoplasmic side. Cytoplasm, cytoskeleton, cilium axoneme. Cytoplasm, cytoskeleton, cilium basal body. Note=Appears in a pattern of punctae flanking the microtubule axoneme that likely correspond to small membrane-associated patches. Localizes to the so-called ciliary gate where vesicles carrying ciliary cargo fuse with the membrane.

DISEASE:

Defects in ARL6 are a cause of Bardet-Biedl syndrome type 3 (BBS3) [MIM:209900]. Bardet-Biedl syndrome (BBS) is a genetically heterogeneous disorder characterized by usually severe pigmentary retinopathy, early onset obesity, polydactyly, hypogenitalism, renal malformation and mental retardation. Secondary features include diabetes mellitus, hypertension and congenital heart disease.

Defects in ARL6 are the cause of retinitis pigmentosa type 55 (RP55) [MIM:613575]. RP55 is a retinal dystrophy belonging to the group of pigmentary retinopathies. Retinitis pigmentosa is characterized by retinal pigment deposits visible on fundus examination and primary loss of rod photoreceptor cells followed by secondary loss of cone photoreceptors. Patients typically have night vision blindness and loss of midperipheral visual field. As their condition progresses, they lose their far peripheral visual field and eventually central vision as well.

Similarity:

Belongs to the small GTPase superfamily. Arf family.

SWISS:

Q9H0F7

Gene ID:

84100

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.