

分离蛋白 SCRN1 抗体

产品货号： mIR19600

英文名称： SCRN1

中文名称： 分离蛋白 SCRN1 抗体

别名： KIAA0193; OTTHUMP00000202313; SCRN 1; SCRN1; SCRN1_HUMAN; Secernin 1; Secernin-1; Secernin1; SES1.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500

(石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 46kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from human SCRN1:151-250/414

亚型： IgG

纯化方法： affinity purified by Protein A

储存液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed

产品介绍： This gene likely encodes a member of the secernin family of proteins. A similar protein in rat functions in regulation of exocytosis in mast cells. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2009]

Function:

Regulates exocytosis in mast cells. Increases both the extent of secretion and the sensitivity of mast cells to stimulation with calcium.

Subcellular Location:

Cytoplasm.

Similarity:

Belongs to the peptidase C69 family. Secernin subfamily.

SWISS:



Q12765

Gene ID:

9805

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.