

## 1 号染色体开放阅读框 31 抗体

产品货号： mlR9785

英文名称： C1orf31

中文名称： 1 号染色体开放阅读框 31 抗体

别名： C1orf31; CA031\_HUMAN; Chromosome 1 open reading frame 31; Hypothetical protein LOC388753; RP5-827C21.3; Uncharacterized protein C1orf31.

研究领域： 细胞生物 免疫学 干细胞

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Pig, Cow, Horse, Sheep,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:50-200 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量：14kDa

细胞定位：细胞浆

性状：Lyophilized or Liquid

浓度：1mg/ml

免疫原：KLH conjugated synthetic peptide derived from human C1orf31:51-125/125

亚型：IgG

纯化方法：affinity purified by Protein A

储存液：0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**保存条件：** Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed：** PubMed

**产品介绍 background：**

The cytochrome c oxidase (COX) family of proteins function as the final electron donor in the respiratory chain to drive a proton gradient across the inner mitochondrial membrane, ultimately resulting in the production of

water. C1orf31 is a 125 amino acid mitochondrial protein that belongs to the cytochrome c oxidase subunit 6B family. There are three isoforms of C1orf31 that are produced as a result of alternative splicing events. The gene encoding C1orf31 maps to human chromosome 1, the largest human chromosome spanning about 260 million base pairs and making up 8% of the human genome. There are about 3,000 genes on chromosome 1, and considering the great number of genes there are also a large number of diseases associated with chromosome 1. Notably, the rare aging disease Hutchinson-Gilford progeria is associated with the LMNA gene which encodes lamin A. When defective, the LMNA gene product can build up in the nucleus and cause characteristic nuclear blebs. The mechanism of rapidly enhanced aging is unclear and is a topic of continuing exploration.

**Subunit:**

Interacts with COA1.

**Subcellular Location:**

Mitochondrion.

**Similarity:**

Belongs to the cytochrome c oxidase subunit 6B family.

**SWISS:**

Q5JTJ3

**Gene ID:**

388753

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic



applications.