

## 磷脂酰肌醇激酶催化亚单位 D 抗体

产品货号： mIR10656

英文名称： PI3 Kinase p110 delta

中文名称： 磷脂酰肌醇激酶催化亚单位 D 抗体

别名： p110D; p110delta; Phosphatidylinositol 3 kinase catalytic delta polypeptide; Phosphatidylinositol 4 5 bisphosphate 3 kinase catalytic subunit delta isoform; Phosphatidylinositol 45 bisphosphate 3 kinase catalytic subunit delta isoform; Phosphoinositide 3 kinase B; Phosphoinositide 3 kinase C; Phosphoinositide 3 kinase catalytic delta polypeptide; PI3 kinase p110 subunit delta; PI3K; PIK3CD; PK3CD; PtdIns 3 kinase p110; PI 3 Kinase p110 delta; PK3CD\_HUMAN; Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit delta isoform; PI3-kinase subunit delta; PI3K-delta; PI3Kdelta; PtdIns-3-kinase subunit delta; Phosphatidylinositol 4,5-bisphosphate 3-kinase 110 kDa catalytic subunit delta; PtdIns-3-kinase subunit p110-delta; p110delta.

研究领域： 肿瘤 免疫学 信号转导 转录调节因子 激酶和磷酸酶

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500

(石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 115kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human PI3 Kinase p110 delta:51-150/1044

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed : PubMed

产品图片

