

## 含 SH2 结构域接头蛋白 B 抗体

产品货号： mIR19757

英文名称： SHB

中文名称： 含 SH2 结构域接头蛋白 B 抗体

别名： bA3J10.2; OTTHUMP00000021397; RP11 3J10.8; RP113J10.8; SH2 domain containing adapter protein B; SH2 domain-containing adapter protein B; SHB (Src homology 2 domain containing) adaptor protein B; SHB adaptor protein (a Src homology 2 protein); SHB adaptor protein; Shb; SHB\_HUMAN; Src homology 2 domain containing adaptor protein B; Src homology 2 protein.

研究领域： 细胞生物 发育生物学 信号转导 干细胞

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 55kDa

细胞定位： 细胞浆 细胞膜

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human SHB:201-300/509

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed :** PubMed