

## 含 SH2 结构域蛋白 D4A 抗体

产品货号： mlR19745

英文名称： SH2D4A

中文名称： 含 SH2 结构域蛋白 D4A 抗体

别 名： FLJ20967; PPP1R38; Protein phosphatase 1 regulatory subunit 38; Protein SH; Protein SH(2)A; Protein SH2A; SH2 domain containing 4A; SH2 domain containing protein 4A; SH2 domain-containing protein 4A; SH24A\_HUMAN; SH2A; Sh2d4a.

研究领域： 肿瘤 细胞生物 信号转导 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Cow, Horse,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量： 53kDa

细胞定位： 细胞浆

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

**免 疫 原：** KLH conjugated synthetic peptide derived from human SH2D4A:201-300/454

**亚 型：** IgG

**纯化方法：** affinity purified by Protein A

**储 存 液：** 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**保存条件：** Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed：** PubMed