

含普列克底物蛋白家族 B 成员 2 抗体

产品货号： mIR12723

英文名称： PLEKHB2

中文名称： 含普列克底物蛋白家族 B 成员 2 抗体

别名： Evectin 2; Evectin-2; Evectin2; EVT2; FLJ20783; FLJ23679; FLJ30436; PH domain containing family B member 2; PH domain-containing family B member 2; PKHB2_HUMAN; Pleckstrin homology domain containing family B member 2; Pleckstrin homology domain containing, family B (evectins) member 2; Pleckstrin homology domain-containing family B member 2; Plekhhb2.

研究领域： 细胞生物 免疫学 信号转导

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Pig, Cow, Horse, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 25kDa

细胞定位： 细胞膜

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human PLEKHB2:131-222/222

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed