

地中海热蛋白 MEFV 抗体

产品货号： mIR18750

英文名称： MEFV

中文名称： 地中海热蛋白 MEFV 抗体

别名： FMF; Marenostrin; Mediterranean fever; Mediterranean fever protein; MEF; Mefv; MEFV_HUMAN; Pypin; TRIM20.

研究领域： 细胞生物 信号转导 结合蛋白 细胞骨架

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 86kDa

细胞定位： 细胞核 细胞浆

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原 : KLH conjugated synthetic peptide derived from human MEFV:501-600/781

亚型 : IgG

纯化方法 : affinity purified by Protein A

储存液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : This gene encodes a protein, also known as pyrin or marenastrin, that is an important modulator of innate immunity. Mutations in this gene are associated with Mediterranean fever, a hereditary periodic fever syndrome. [provided by RefSeq, Jul 2008]

Function:

Involved in innate immunity and the inflammatory response. Interacts with several components of the inflammasome complex, a large oligomeric structure which recruits and activates CASP1 and ultimately induces maturation of cytokines such as IL1B. However, the exact role of MEFV in the inflammatory pathway is uncertain as contradictory effects on IL1B processing have been reported in different experimental systems. Has been shown to activate IL1B production (Ref.13). Has also been shown to inhibit IL1B production (Ref.14, Ref.15). Also required for PSTPIP1-induced PYCARD oligomerization and for formation of pyroptosomes, large supramolecular structures composed of oligomerized PYCARD dimers which form prior to inflammatory apoptosis. Can reduce PYCARD-induced apoptosis. Recruits PSTPIP1 to pyroptosomes, and required for PSTPIP1 oligomerization.

Subcellular Location:

Nucleus and Cytoplasm > cytoskeleton. Associated with microtubules and with the filamentous actin of perinuclear filaments and peripheral lamellar ruffles.

Tissue Specificity:

Expressed in peripheral blood leukocytes, particularly in mature granulocytes and to a lesser extent in monocytes but not in lymphocytes. Detected in spleen, lung and muscle, probably as a result of leukocyte infiltration in these tissues. Not expressed in thymus, prostate, testis, ovary, small intestine, colon, heart, brain, placenta, liver, kidney, pancreas. Expression detected in several myeloid leukemic, colon cancer, and prostate cancer cell lines.

Post-translational modifications:

Cleaved by CASP1 Probable. The N-terminal cleavage product localizes to the nucleus as a filamentous network and to the cytoplasm, interacts more strongly with RELA and NFKBIA than the full-length protein, enhances the nuclear localization of RELA and induces NFKBIA proteolysis. The C-terminal cleavage product localizes to the cytoplasm.

DISEASE:

Defects in MEFV are the cause of familial Mediterranean fever autosomal recessive (ARFMF) [MIM:249100]. ARFMF is an inherited disorder characterized by recurrent episodic fever, serosal inflammation and pain in the abdomen, chest or joints. ARFMF is frequently complicated by amyloidosis, which leads to renal failure and can be prophylactically treated with colchicine. ARFMF primarily affects ancestral ethnic groups living around the Mediterranean basin: North African Jews, Armenians, Arabs and Turks. The disease is also distributed in other populations including Greeks, Cypriots, Italians and Spanish, although at a lower prevalence.

Defects in MEFV are the cause of familial Mediterranean fever autosomal dominant (ADFMF) [MIM:134610]. ADFMF is characterized by periodic fever, serosal inflammation and pain in the abdomen, chest or joints as seen also in the autosomal recessive form of the disease. It is associated with renal amyloidosis and characterized by colchicine unresponsiveness.

Similarity:

Contains 1 B box-type zinc finger.

Contains 1 B30.2/SPRY domain.

Contains 1 DAPIN domain.



SWISS:

O15553

Gene ID:

4210

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.