

DNA 聚合酶 α 2 抗体

产品货号 : mlR12754

英文名称 : POLA2

中文名称 : DNA 聚合酶 α 2 抗体

别 名 : AI573378; DNA directed DNA polymerase alpha 2; DNA polymerase alpha 70 kDa subunit; DNA polymerase alpha subunit B; DNA polymerase alpha subunit II; DNA polymerase subunit alpha B; DNA polymerase subunit II; DPOA2_HUMAN; FLJ21662; FLJ37250; polA2; Polymerase (DNA directed) alpha 2 (70kD subunit); Polymerase (DNA directed), alpha (70kD); Polymerase (DNA directed), alpha 2 accessory subunit.

研究领域 : 细胞生物 细胞周期蛋白 表观遗传学

抗体来源 : Rabbit

克隆类型 : Polyclonal

交叉反应 : Human, Mouse, Rat, Dog, Pig, Horse, Rabbit, Zebrafish,

产品应用 : IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 (石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 66kDa

细胞定位 : 细胞核

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from human POLA2:101-200/598

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed