

多腺苷二磷酸多聚酶 PARP8 抗体

产品货号： mlR19889

英文名称： PARP8

中文名称： 多腺苷二磷酸多聚酶 PARP8 抗体

别名： ADP-ribosyltransferase diphtheria toxin-like 14; ARTD14; AW558171; DDF1; DKFZp434J214; DKFZp686N0351; DKFZp686P1838; FLJ40466; PARP 1; PARP 7; PARP-7; PARPT_HUMAN; pART14; Poly [ADP ribose] polymerase 7; Poly [ADP-ribose] polymerase 7; TCDD inducible poly(ADP ribose) polymerase; TCDD-inducible poly [ADP-ribose] polymerase; TIPARP.

研究领域： 细胞生物 结合蛋白 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Pig, Cow, Horse, Rabbit, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 96kDa

细胞定位： 细胞核 细胞浆

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human PARP8:701-800/854

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed