

线粒体核糖体蛋白 L46 抗体

产品货号： mIR17797

英文名称： MRPL46

中文名称： 线粒体核糖体蛋白 L46 抗体

别名： 39S ribosomal protein L46; 39S ribosomal protein L46 mitochondrial; C15orf4; L46mt; LIECG 2; LIECG2; MGC22762; mitochondrial; Mitochondrial ribosomal protein L46; MRP L46; MRP-L46; MRPL 46; Mrpl46; P2ECSL; RM46_HUMAN.

研究领域： 细胞生物 转录调节因子 线粒体 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Pig, Cow, Horse, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需

做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 32kDa

细胞定位 : 细胞浆 线粒体

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human MRPL46:31-130/279

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. [provided by RefSeq, Jul 2008]

Subcellular Location:

Mitochondrion.

Similarity:

Belongs to the ribosomal protein L46 family.

SWISS:

Q9H2W6

Gene ID:

26589

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

