

膜蛋白 CLDND2 抗体

产品货号： mlR8014

英文名称： CLDND2

中文名称： 膜蛋白 CLDND2 抗体

别名： Claudin domain containing 2; Claudin domain containing protein 2; PMP22CD.

研究领域： 细胞生物 免疫学 细胞粘附分子

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Pig, Cow, Sheep,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量：18kDa

细胞定位：细胞膜

性状：Lyophilized or Liquid

浓度：1mg/ml

免疫原：KLH conjugated synthetic peptide derived from human CLDND2:101-167/167

亚型：IgG

纯化方法：affinity purified by Protein A

储存液：0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件：Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed：PubMed

产品介绍：CLDND2 belongs to the PMP-22/EMP/MP20 family of proteins. At present the molecular function is unknown.

Subcellular Location:

Membrane; Multi-pass membrane protein (Potential).

Similarity:

Belongs to the PMP-22/EMP/MP20 family.

SWISS:

Q8NHS1

Gene ID:

125875

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品图片

