

20 号染色体开放阅读框 50 抗体

产品货号： mlR4457

英文名称： C20orf50

中文名称： 20 号染色体开放阅读框 50 抗体

别名： SEL1L2; C20orf50; Protein sel-1 homolog 2; SE1L2_HUMAN; Sel 1 suppressor of lin 12 like 2 (C. elegans); Sel-1L2; Sel1l2; Suppressor of lin-12-like protein 2.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Rabbit,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500
(石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量：75kDa

细胞定位：细胞膜

性状：Lyophilized or Liquid

浓度：1mg/ml

免疫原：KLH conjugated synthetic peptide derived from human C20orf50:131-230/688

亚型：IgG

纯化方法：affinity purified by Protein A

储存液：Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4

保存条件：Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed：PubMed

产品介绍 background:

SEL1L2 (SEL1L2 ERAD E3 Ligase Adaptor Subunit) is a Protein Coding gene. Among its related pathways are Protein processing in endoplasmic reticulum. An important paralog of this gene is SEL1L.

Subcellular Location:

Membrane.

Similarity:

Belongs to the sel-1 family.

Contains 11 Sel1-like repeats.

SWISS:

Q5TEA6

Gene ID:

80343

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品图片

