

Rab27a 蛋白抗体

产品货号： mlR1581

英文名称： Rab27a

中文名称： Rab27a 蛋白抗体

别名： GS2; GTP-binding protein Ram; HsT18676; MGC117246; Rab-27; RAB-27A; RAB27; RAB27A; RAB27A member RAS oncogene family; RAM; Ras-related protein Rab-27A; Ras-related protein Rab27A; RB27A_HUMAN.

研究领域： 肿瘤 细胞生物 转录调节因子 合成与降解 转运蛋白

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Pig, Cow, Horse, Sheep,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 25kDa

细胞定位： 细胞浆 细胞膜

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from human RAM-11:179-217/221

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : Summary: The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. The protein is membrane-bound and may be involved in protein transport and small GTPase mediated signal transduction. Mutations in this gene are associated with Griscelli syndrome type 2. Alternative splicing occurs at this locus and four transcript variants encoding the same protein have been identified. [provided by RefSeq].

Function:

Plays a role in cytotoxic granule exocytosis in lymphocytes. Required for both granule maturation and granule docking and priming at the immunologic synapse.

Subunit:

Binds SYTL1, SLAC2B, MYRIP, SYTL3, SYTL4 and SYTL5 (By similarity). Binds MLPH and SYTL2. Interacts with UNC13D.

Subcellular Location:

Membrane; Lipid-anchor. Melanosome. Late endosome. Lysosome. Note=Identified by mass spectrometry in melanosome fractions from stage I to stage IV. Localizes to endosomal exocytic vesicles.

Tissue Specificity:

Found in all the examined tissues except in brain. Low expression was found in thymus, kidney, muscle and placenta. Detected in melanocytes, and in most tumor cell lines examined. Expressed in cytotoxic T-lymphocytes (CTL) and mast cells.

DISEASE:

Griscelli syndrome 2 (GS2) [MIM:607624]: Rare autosomal recessive disorder that results in pigmentary dilution of the skin and hair, the presence of large clumps of pigment in hair shafts, and an accumulation of melanosomes in melanocytes. GS2 patients also develop an uncontrolled T-lymphocyte and macrophage activation syndrome, known as hemophagocytic syndrome, leading to death in the absence of bone marrow transplantation. Neurological impairment is present in some patients, likely as a result of hemophagocytic syndrome. Note=The disease is caused by mutations affecting the gene represented in this entry.

Similarity:

Belongs to the small GTPase superfamily. Rab family.

SWISS:

P51159

Gene ID:

5873

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.