

TOMM20L 蛋白抗体

产品货号： mlR17224

英文名称： TOMM20L

中文名称： TOMM20L 蛋白抗体

别名： 4930553D19Rik; TO20L_HUMAN; TOMM20-like protein 1; Tomm20l; Translocase of outer mitochondrial membrane 20 homolog (yeast) like.

研究领域： 细胞生物 转录调节因子 细胞类型标志物 线粒体

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 18kDa

细胞定位 : 细胞浆 线粒体

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human TOMM20L:21-120/152

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.



PubMed : PubMed