

## 鞘脂激活蛋白样蛋白 1 抗体

产品货号： mIR19450

英文名称： PSAPL1

中文名称： 鞘脂激活蛋白样蛋白 1 抗体

别 名： Proactivator polypeptide like 1; Prosaposin like 1; Prosaposin like protein 1; PSAPL1; SAPL1\_HUMAN; Saposin A like; Saposin B like; Saposin B Val like; Saposin C like; Saposin D like; Saposin D-like.

研究领域： 肿瘤 细胞生物 神经生物学 信号转导

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat,

产品应用： IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量： 57kDa

细胞定位： 分泌型蛋白

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from human PSAPL1:51-150/521

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed :** PubMed