

## 地高辛抗体

产品货号： mIR0356

英文名称： Digoxin

中文名称： 地高辛抗体

别 名： Digoxigenin; Lanoxin.

产品类型： 药物与化合物抗体

研究领域： 细胞生物 免疫学 药物及化合物

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Digoxin

产品应用： ELISA=1:500-1000

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量： 0.78kDa

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

免 疫 原： BSA conjugated to Dig:

亚 型： IgG

纯化方法： affinity purified by Protein A

**储 存 液：** 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**保存条件：** Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed：** PubMed

**产品介绍：** Digoxin is a beta blocker drug originally derived from the foxglove plant, Digitalis lanata. Digoxin is used primarily to improve the pumping ability of the heart in congestive heart failure (CHF), and treat problems such as high blood pressure. It is also used to help normalize some dysrhythmias (abnormal types of heartbeat).

**SWISS:**

N/A

**CAS:**

20830-75-5

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.