

核输入蛋白 NIP7 抗体

产品货号： mlR19249

英文名称： NIP7

中文名称： 核输入蛋白 NIP7 抗体

别名： 60S ribosome subunit biogenesis protein NIP7 homolog; CGI 37; CGI37; FLJ10296; HSPC031; HSPC180; KD93; NIP 7; NIP7; NIP7 nucleolar pre rRNA processing protein; NIP7_HUMAN; Nuclear import 7; Nuclear import 7 homolog (S. cerevisiae); Nuclear import 7 homolog; OK/SW cl.76; OK/SW cl.78; OTTHUMP00000174889.

研究领域： 染色质和核信号 干细胞 转录调节因子 结合蛋白 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Chicken, Dog, Pig, Cow, Horse, Rabbit, Zebrafish, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 20kDa

细胞定位： 细胞核

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human NIP7:61-160/180

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed