

## 嗅觉受体 2A1 抗体

产品货号： mIR19708

英文名称： OR2A1

中文名称： 嗅觉受体 2A1 抗体

别名： Olfactory receptor 2A1/2A42; Olfactory receptor family 2 subfamily A member 1; Olfactory receptor family 2 subfamily A member 42; Olfactory receptor OR7-16; Olfactory receptor OR7-19; OR2A1; OR2A1\_HUMAN; OR2A42.

研究领域： 细胞生物 神经生物学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 35kDa

细胞定位： 细胞膜

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human OR2A1:221-310/310 <Cytoplasmic>

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed :** PubMed