

SAR1 蛋白抗体

产品货号： mlR19391

英文名称： SAR1A

中文名称： SAR1 蛋白抗体

别名： COPII associated small GTPase; COPII-associated small GTPase; GTP binding protein SAR1a; GTP-binding protein SAR1a; Masra 2; Masra2; SAR 1; SAR1 gene homolog A; SAR1 homolog A (S. cerevisiae); SAR1; SAR1a gene homolog (S cerevisiae) 1; SAR1a gene homolog 1 (S cerevisiae); SAR1a gene homolog 1; SAR1A_HUMAN; SARA 1; Sara.

研究领域： 信号转导 结合蛋白 G 蛋白信号

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Pig, Cow, Xenopus laevis

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 22kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human SAR1:51-150/198

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed