

热休克蛋白 70 蛋白 2 抗体

产品货号： mlR18080

英文名称： HSPA2

中文名称： 热休克蛋白 70 蛋白 2 抗体

别名： 70kDa; DAQB 147D11.1 001; Hcp70.2; Heat shock 70 kDa protein 2; heat shock 70kDa protein 1A; Heat shock protein 2; heat shock protein 70; Heat shock protein 70.2; Heat shock related 70 kDa protein 2; Heat shock-related 70 kDa protein 2; Heat-shock protein, 70-KD, 2; Heat-shock protein, 70-KD, 3; HSP70 2; HSP70 3; Hsp70-2; HSP70-3; HSP70.2; HSP70A2; HSP72; HSP72_HUMAN; HSPA2; Hspt70; Hst70; MGC58299; MGC7795; MGC93458; OTTHUMP00000180664; Testis-specific heat shock protein-related; XXbac BCX40G17.3 001.

研究领域： 肿瘤 细胞生物 信号转导

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Pig, Cow, Horse, Rabbit, Sheep, Monkey,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500
(石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 70kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human HSPA2:451-550/639

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed