

## SAMD4B 蛋白抗体

产品货号： mlR19375

英文名称： SAMD4B

中文名称： SAMD4B 蛋白抗体

别名： FLJ10211; HSmaug2; MGC99832; Protein Smaug homolog 2; SAM domain-containing protein 4B; Smaug 2; Smaug homolog B; SMAUG2; SMGB; Sterile alpha motif domain containing 4B; Sterile alpha motif domain-containing protein 4B.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Pig, Cow, Horse, Rabbit, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 75kDa

细胞定位： 细胞核 细胞浆

性状： Lyophilized or Liquid

浓度： 1mg/ml

**免 疫 原：** KLH conjugated synthetic peptide derived from human SAMD4B:231-330/694

**亚 型：** IgG

**纯化方法：** affinity purified by Protein A

**储 存 液：** 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**保存条件：** Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed：** PubMed

**产品介绍：** Sterile alpha motif is a widespread domain in signalling and nuclear proteins. In EPH-related tyrosine kinases, it appears to mediate cell-cell initiated signal transduction via the binding of SH2-containing proteins to a conserved tyrosine. SAMD4B (Sterile alpha motif domain-containing protein 4B) belongs to the SMAUG family. SAMD4B contains 1 SAM (sterile alpha motif) domain.

**Subcellular Location:**

Cytoplasmic and Nuclear

**SWISS:**

Q5PRF9

**Gene ID:**

55095

**Important Note:**



This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.