

热休克因子结合蛋白 1 样蛋白 1 抗体

产品货号： mIR17397

英文名称： HSBP1L1

中文名称： 热休克因子结合蛋白 1 样蛋白 1 抗体

别名： FLJ10967; Heat shock factor binding protein 1 like 1; MGC189743.

研究领域： 细胞生物 转录调节因子

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Rat, Bee,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 8kDa

细胞定位 : 细胞核

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human HSBP1L1:131-230/342

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed