

子宫珠蛋白抗体

产品货号： mlR1935

英文名称： Uteroglobin

中文名称： 子宫珠蛋白抗体

别名： blastokinin; CC10; CC16; CCPBP; CCSP; Clara cell phospholipid binding protein; Clara cell phospholipid-binding protein; clara cell secretory protein; Clara cell specific 10 kD protein; Clara Cells 10 kDa Secretory Protein; OTTHUMP00000236107; SCGB1A1; Secretoglobin family 1A member 1; Secretoglobin, family 1A, member 1 (uteroglobin); UG; UGB; UP-1; UP1; Urinary Protein 1; Urine Protein 1; UTER_HUMAN; Uteroglobin; Utg.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Mouse, Rat,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 8kDa

细胞定位： 分泌型蛋白

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原 : KLH conjugated synthetic peptide derived from mouse Uteroglobin:51-96/96

亚型 : IgG

纯化方法 : affinity purified by Protein A

储存液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : Uteroglobin belongs to the family of secretoglobins and is a secreted protein product of nonciliated bronchiolar Clara cells. There is convincing data suggesting it has phospholipase A2 inhibitory activity, as well as, a number of other immunomodulatory features including inhibition of interferon gamma signalling and Th1 vs. Th2 lymphocyte regulation. It was proposed as a potential peripheral marker of respiratory epithelial injury and bronchial dysfunction. Clara Cell Protein 16 concentrations have been determined in both serum and bronchoalveolar lavage fluid in numerous studies since 1994. In serum, its increase is associated with age, asbestos, nitrogen chloride and ozone exposure, sarcoidosis and high PEEP ventilation. Decreased serum CC16 levels are found after pulmonary resection, in silica-exposed workers, smokers and in asthma. Decreased CC16 concentrations were also found in the amniotic fluid of fetuses suffering from pulmonary hypoplasia caused by various mechanisms (diaphragmatic hernia, diabetic fetopathy, Turner and Down syndrome). In pleural effusions, the CC16 concentration appears to be associated with its diffusion from the lung as evidenced by high CC16 levels in cardiac pleural congestion.

Function:

Binds phosphatidylcholine, phosphatidylinositol, polychlorinated biphenyls (PCB) and weakly progesterone, potent inhibitor of phospholipase A2.

Subunit:

Homodimer; antiparallel disulfide-linked.

Subcellular Location:

Secreted.

Tissue Specificity:

Clara cells (nonciliated cells of the surface epithelium of the pulmonary airways).

Similarity:

Belongs to the secretoglobin family.

SWISS:

Q06318

Gene ID:

22287

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

子宫珠蛋白(uteroglobin)是一种类固醇诱导的类似细胞因子的多功能蛋白质，具有抗炎和免疫调节特性,它不仅在子宫内表达,而且可以在甲状腺,呼吸道,消化道，乳腺、胰腺、前列腺、脑垂体、睾丸等都有不同的表达。

mlbio 酶联生物
Good elisakit producers