

艾滋病病毒 EP2 抗体

产品货号： mIR17460

英文名称： HIVEP2

中文名称： 艾滋病病毒 EP2 抗体

别名： HIV EP2; HIV-EP2; HIVEP2; human immunodeficiency virus type I enhancer binding protein 2; Human immunodeficiency virus type I enhancer-binding protein 2; MBP 2; MBP-2; MHC binding protein 2; MHC-binding protein 2; MIBP1; SHN2; Transcription factor HIVEP2; ZAS2; ZEP2_HUMAN; ZNF40B.

研究领域： 细胞生物 免疫学 淋巴细胞 t-淋巴细胞 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Chicken, Dog, Pig, Cow, Horse, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需

做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 269, 35kDa

细胞定位 : 细胞核

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human HIVEP2:321-420/2446

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.



PubMed : PubMed