

锌指蛋白 202 抗体

产品货号： mlR12799

英文名称： ZNF202

中文名称： 锌指蛋白 202 抗体

别名： Zinc finger protein 202; Zinc finger protein with KRAB and SCAN domains 10; ZKSCAN 10; ZKSCAN10; ZN202_HUMAN; ZNF 202; ZNF202.

研究领域： 细胞生物 转录调节因子 锌指蛋白 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Pig, Cow, Horse, Rabbit, Sheep,

产品应用： IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）
not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

分子量： 75kDa

细胞定位： 细胞核

性状： Lyophilized or Liquid

浓度： 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human ZNF202:1-100/648

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed