

## 锌指蛋白 215 抗体

产品货号： mlR12802

英文名称： ZNF215

中文名称： 锌指蛋白 215 抗体

别名： BAZ 2; BAZ-2; BAZ2; BWSCR2 associated zinc finger protein 2; BWSCR2-associated zinc finger protein 2; Zinc finger protein 215; Zinc finger protein with KRAB and SCAN domains 11; ZKSCAN 11; ZKSCAN11; ZN215\_HUMAN; ZNF 215; ZNF215.

研究领域： 转录调节因子 锌指蛋白 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 60kDa

细胞定位： 细胞核

性状： Lyophilized or Liquid

浓度： 1mg/ml

**免 疫 原：** KLH conjugated synthetic peptide derived from human ZNF215:301-400/517

**亚 型：** IgG

**纯化方法：** affinity purified by Protein A

**储 存 液：** 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**保存条件：** Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed：** PubMed