

锌指蛋白 248 抗体

产品货号： mlR12804

英文名称： ZNF248

中文名称： 锌指蛋白 248 抗体

别名： bA162G10.3 (zinc finger protein); bA162G10.3; KIAA0137; KRAB protein domain; Kruppel-type zinc finger protein; MGC131882; OTTHUMP00000019462; OTTHUMP00000019463; OTTHUMP00000019464; Zinc finger protein 248; ZN248_HUMAN; ZNF248.

研究领域： 转录调节因子 锌指蛋白 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 67kDa

细胞定位： 细胞核

性状： Lyophilized or Liquid

浓度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from human ZNF248:131-230/579

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed