

锌指蛋白 313 抗体

产品货号： mlR12806

英文名称： ZNF313

中文名称： 锌指蛋白 313 抗体

别 名： 1110008J21Rik; AI225886; AW549494; MGC112560; OTTMUSP00000001082;
OTTMUSP00000017392; PSORS12; RING finger protein 114; RN114_HUMAN; RNF114; RP23-118A2.1; Zfp228;
Zfp313; Zinc finger protein 228; Zinc finger protein 313; ZNF228; ZNF313。

研究领域： 细胞生物 转录调节因子 锌指蛋白 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Chicken, Dog, Pig, Cow, Horse, Rabbit, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需
做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量： 26kDa

细胞定位： 细胞核 细胞浆

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from human ZNF313:101-200/228

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed