

# 肌肉特定环指蛋白 3 抗体

产品货号： mlR9376

英文名称： MURF3

中文名称： 肌肉特定环指蛋白 3 抗体

别名： MURF 3; MURF; MURF-3; MuRF3; Muscle specific RING finger protein 3; Muscle specific RING finger protein; Muscle specific RING finger protein homolog; Muscle-specific RING finger protein 3; Muscle-specific RING finger protein; Ring finger protein 30; RNF30; TRI54\_HUMAN; TRIM 54; TRIM54; RNF28; Tripartite motif containing 54; tibody Tripartite motif-containing protein 54.

研究领域： 细胞生物 免疫学 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Pig, Cow, Rabbit,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:50-200 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 40kDa

细胞定位 : 细胞浆

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human MURF3:131-240/358

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed :** PubMed

**产品介绍：** background

MuRF1, is a nuclear protein that interacts with SMT3b and the large myofibrillar protein Titin. In muscle cells, MuRF2 (RNF29) regulates gene expression and protein turnover. It localizes to the cytoplasm, but under atrophic conditions it is detected in the nucleus. MuRF2 can form oligomers with various other proteins, including Titin and Myosin. MuRF3, also designated tripartite motif-containing 54 (TRIM54) or ring finger protein 30 (RNF30), interacts with tubulin and stabilizes microtubules during myotube formation. It is a cytoplasmic protein that localizes to the Z-lines in skeletal muscles, while MuRF2 localizes to the sarcomeric M-band in cardiomyocytes. MuRF3 shares 77% and 65% sequence identity with MuRF1 and MuRF2, respectively. MuRF1-3 share a conserved N-terminal RING domain and zinc-binding B-box motif, and two coiled-coil dimerization motif boxes, in their central regions.

**Function:**

May bind and stabilize microtubules during myotubes formation (By similarity).

**Subunit:**

Homooligomer and heterooligomer. Interacts with tubulin. Interacts with TRIM63 and probably with TRIM55.

**Subcellular Location:**

Cytoplasm

**Tissue Specificity:**

Specifically expressed in heart and skeletal muscle.

**Similarity:**

Contains 1 B box-type zinc finger.

Contains 1 COS domain.

Contains 1 RING-type zinc finger.

**SWISS:**

Q9BYV2

**Gene ID:**

57159

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.