

鸡白痢、鸡伤寒抗体

产品货号： mlR1675

英文名称： fowl typhoid and pullorum disease

中文名称： 鸡白痢、鸡伤寒抗体

研究领域： 免疫学 细菌及病毒

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Chicken,

产品应用： ELISA=1:500-1000 IF=1:100-500

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

免 疫 原： Inactivated Salmonella gallinarum:

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When

reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : 鸡伤寒和鸡白痢是由鸡白痢沙门氏菌(*salmonella pullorum*)和鸡伤寒沙门氏菌(*salmonella gallinarum*)引起的各种日龄鸡和火鸡的传染病,以急性败血性经过或以慢性隐性感染为特征。

鸡白痢鸡伤寒病依据病原体不同分为三种疾病:鸡白痢、鸡伤寒和副伤寒。临床上多表现为败血症和肠炎,也是禽类最为重要的细菌蛋传性疾病之一,导致种禽场和商品禽场的生产性能下降。