

鸡传染性喉气管炎病毒糖蛋白 E 抗体

产品货号： mlR1674

英文名称： ILTV glycoprotein E

中文名称： 鸡传染性喉气管炎病毒糖蛋白 E 抗体

别名： Infections Laryngotracheitis Virus; ILTV glycoprotein E; Infections Laryngotracheitis Virus glycoprotein E; US8.

研究领域： 细菌及病毒

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： ILTV

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 55kDa

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from Infections Laryngotracheitis Virus glycoprotein E:41-150/499

亚型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed

产品介绍： 鸡传染性喉支气管炎是由疱疹病毒科的传染性喉支气管炎病毒引起的一种急性呼吸道传染病，主要以呼吸困难，咳嗽，咳出含有血液的渗出物，喉部和气管黏膜肿胀、出血，糜烂等症状为特征。