

## 甲基 CpG 结合域蛋白 3 样 3 抗体

产品货号： mIR18703

英文名称： MBD3L3

中文名称： 甲基 CpG 结合域蛋白 3 样 3 抗体

别名： MB3L3\_HUMAN; MBD3 like 3; MBD3 like protein 3; MBD3-like protein 3; MBD3L3; putative methyl CpG binding domain protein 3 like 3; Putative methyl-CpG-binding domain protein 3-like 3.

研究领域： 转录调节因子 结合蛋白 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 23kDa

细胞定位： 细胞核

性状： Lyophilized or Liquid

浓度： 1mg/ml

**免疫原** : KLH conjugated synthetic peptide derived from human MBD3L3:81-180/208

**亚型** : IgG

**纯化方法** : affinity purified by Protein A

**储存液** : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**保存条件** : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed** : PubMed