

N- α -乙酰转移酶 40 抗体

产品货号： mIR18990

英文名称： NAA40

中文名称： N- α -乙酰转移酶 40 抗体

别名： N acetyltransferase 11 (GCN5 related, putative); N alpha acetyl transferase 40; N alpha acetyltransferase 40; N(alpha) acetyltransferase 40; N(alpha) acetyltransferase 40, NatD catalytic subunit, homolog; N(alpha)-acetyltransferase 40, NatD catalytic subunit, homolog (S. cerevisiae); N-acetyltransferase 11; N-alpha-acetyltransferase 40, NatD catalytic subunit; naa40; NAA40_HUMAN; NAT11; NatD catalytic subunit; PATT1.

研究领域： 肿瘤 细胞生物 细胞凋亡

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Pig, Cow, Horse, Sheep,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 24kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human NAA40:51-150/237

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed