

胰岛素样肽 3 抗体

产品货号： mIR16606

英文名称： DILP3

中文名称： 胰岛素样肽 3 抗体

别 名： CG14167; DILP 3; DILP; DILP3; DmelCG14167; Ilp3; INSL3_DROME; Insulin like peptide 3; Insulin-related peptide 3; lip; Probable insulin-like peptide 3 B chain.

研究领域： 细胞生物 信号转导 生长因子和激素

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Mouse, Rat, Firefly, Drosophila melanogaster

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量： 14kDa

细胞定位： 分泌型蛋白

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from Fire fly DILP3:21-100/120

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed : PubMed

产品图片

