

荆豆凝集素 1 抗体

产品货号： mlR3584

英文名称： Ulex Europaeus Lectin 1

中文名称： 荆豆凝集素 1 抗体

别名： UEA 1; UEA I; UEA-I; UEAI; UEA1; UEA-1; UEA1; Ulex europaeus agglutinin; Anti H(O) lectin 1; Anti H(O) lectin I.

研究领域： 肿瘤 细胞生物 免疫学 细胞凋亡 细胞表面分子 糖蛋白

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Ulexeuropaeus

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:100-500 (石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 27kDa

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from Ulex Europaeus Lectin 1:151-243/243

亚型： IgG

纯化方法： affinity purified by Protein A

储存液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 : Ulex Europaeus is a European gorse shrub with fragrant golden-yellow flowers. Ulex Europaeus-I Lectin is a 46 kDa glycoprotein known to interact with -L fucosyl residues in oligosaccharides present on the membranes of human blood group O erythrocytes, human endothelial cells and a variety of human and animal epithelial cells.

Subunit:

Belongs to the leguminous lectin family.

SWISS:

Q9BTC0

Gene ID:

N/A

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.