

癫痫样相关蛋白 6 抗体

产品货号： mlR6103

英文名称： SEZ6L

中文名称： 癫痫样相关蛋白 6 抗体

别名： Seizure 6 like protein; Seizure related 6 homolog (mouse) like; Seizure related gene 6 (mouse) like; Seizure related gene 6 like; SE6L1_HUMAN.

研究领域： 免疫学 神经生物学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Chicken, Dog, Pig, Cow, Horse, Sheep,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量：109kDa

细胞定位：细胞浆 细胞膜

性状：Lyophilized or Liquid

浓度：1mg/ml

免疫原：KLH conjugated synthetic peptide derived from human SEZ6L:701-800/1024

亚型：IgG

纯化方法：affinity purified by Protein A

储存液：0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件：Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed：PubMed

产品介绍：SEZ6L is a candidate tumor suppressor gene. This protein may contribute to specialized endoplasmic reticulum functions in neurons.

Function:

May contribute to specialized endoplasmic reticulum functions in neurons

Subcellular Location:

Endoplasmic reticulum membrane; Single-pass type I membrane protein

Tissue Specificity:

Widely expressed, including adult and fetal brains and lungs. Not expressed in all lung cancer cell lines.

Similarity:

Belongs to the SEZ6 family.

Contains 3 CUB domains.

Contains 5 Sushi (CCP/SCR) domains.

SWISS:

Q9BYH1

Gene ID:

23544

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品图片

