

## NBPF16 蛋白抗体

产品货号： mIR19034

英文名称： NBPF16

中文名称： NBPF16 蛋白抗体

别名： AG3; Neuroblastoma breakpoint family member 16.

研究领域： 肿瘤 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 78kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from human NBPF16:571-670/670

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed :** PubMed

产品介绍 : Belongs to the NBPF family. Contains 6 NBPF domains.

**Subcellular Location:**

Cytoplasm

**SWISS:**

Q5SXJ2

**Gene ID:**

728936

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.