

MS4A6E 蛋白抗体

产品货号： mIR17843

英文名称： MS4A6E

中文名称： MS4A6E 蛋白抗体

别 名： Membrane-spanning 4-domains subfamily A member 6E; MS4A6E.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量： 16kDa

细胞定位： 细胞膜

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from human MS4A6E:21-100/147 <Extracellular>

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20 ° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed

产品介绍： This gene encodes a member of the membrane-spanning 4A gene family. Members of this nascent protein family are characterized by common structural features and similar intron/exon splice boundaries and display unique expression patterns among hematopoietic cells and nonlymphoid tissues. The gene encoding this protein is localized to 11q12.3, among a cluster of family members. [provided by RefSeq, Jul 2008]

Function:

May be involved in signal transduction as a component of a multimeric receptor complex.

Subcellular Location:

Cell Membrane

Tissue Specificity:

Expressed by malignant and fetal tissue at very low levels.

Similarity:

Similarity: Belongs to the MS4A family

SWISS:

Q96DS6

Gene ID:

245802

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.