

乙型肝炎病毒 X 抗原激活蛋白 5 抗体

产品货号： mIR16959

英文名称： KIAA0020

中文名称： 乙型肝炎病毒 X 抗原激活蛋白 5 抗体

别名： HA-8; HBV X transactivated gene 5 protein; HBV X-transactivated gene 5 protein; HBV XAg-transactivated protein 5; HLA HA8; HLA-HA8; hPUF-A; K0020_HUMAN; Kiaa0020; KIAA0020 protein; MGC8749; Minor histocompatibility antigen HA 8; Minor histocompatibility antigen HA-8; PEN; Penguin homolog; Protein 5 transactivated by hepatitis B virus X antigen (HBxAg); Protein 5 transactivated by hepatitis B virus X antigen; PUF 6; PUF6; Pumilio domain containing protein KIAA0020; Pumilio domain-containing protein KIAA0020; RP11 526D20.2; XTP 5; XTP5.

研究领域： 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Dog, Cow, Rabbit, Sheep, Chimpanzee, Macaque Monkey, Rhesus monkey, Gorilla, Marmoset (common)

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 74kDa

细胞定位： 细胞核

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from human KIAA0020:65-160/648

亚型： IgG

纯化方法： affinity purified by Protein A

储存液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When

reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed