

精子发生相关蛋白 21 抗体

产品货号： mlR17626

英文名称： SPATA21

中文名称： 精子发生相关蛋白 21 抗体

别名： RP11-430L17.3; Spata21; Spergen 2; Spergen2; Spermatogenesis-associated protein 21; SPT21_HUMAN.

研究领域： 细胞生物 发育生物学 干细胞

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Horse,

产品应用： WB=1:500-2000 ELISA=1:500-1000 ICC=1:100-500 IF=1:100-500

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 52kDa

细胞定位： 细胞外基质

性状： Lyophilized or Liquid

浓度： 1mg/ml

免疫原： KLH conjugated synthetic peptide derived from human SPATA21:181-280/469

亚型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

PubMed： PubMed

产品图片

