

## ras 基因相关蛋白 Rab33a 抗体

产品货号： mlR19711

英文名称： RAB33A

中文名称： ras 基因相关蛋白 Rab33a 抗体

别名： MGC1488; RAB 33A; Rab33a; RAB33A member RAS oncogene family; RabS 10; RabS10; Ras related protein Rab 33A; Ras related protein Rab33A; Ras-related protein Rab-33A; RB33A\_HUMAN; Small GTP binding protein S10; Small GTP-binding protein S10.

研究领域： 细胞生物 信号转导 G 蛋白偶联受体 G 蛋白信号

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 27kDa

细胞定位： 细胞膜

性状： Lyophilized or Liquid

浓度： 1mg/ml

**免 疫 原：** KLH conjugated synthetic peptide derived from human RAB33A:1-100/237

**亚 型：** IgG

**纯化方法：** affinity purified by Protein A

**储 存 液：** 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**保存条件：** Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed：** PubMed

**产品介绍：** The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. It is GTP-binding protein and may be involved in vesicle transport. [provided by RefSeq, Jul 2008]

**Subcellular Location:**

Cell membrane.

**Tissue Specificity:**

Expressed only in lymphoid cell lines.

**Similarity:**

Belongs to the small GTPase superfamily. Rab family.

**SWISS:**

Q14088

**Gene ID:**

9363

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.