

溶血磷脂酶样蛋白 1 抗体

产品货号： mlR18580

英文名称： LYPLAL1

中文名称： 溶血磷脂酶样蛋白 1 抗体

别 名： BC027340; FLJ99730; KIAA1238; LYPL1_HUMAN; LYPLAL1; Lysophospholipase like 1; Lysophospholipase-like protein 1; MGC28394; OTTHUMP00000035229.

研究领域： 肿瘤 细胞生物 信号转导 新陈代谢

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Macaque Monkey, Orangutan

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量： 26kDa

细胞定位： 细胞浆

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from human LYPLAL1:61-160/237

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed