

Ras 蛋白脑组织样蛋白 1 抗体

产品货号： mIR19869

英文名称： RHEBL1

中文名称： Ras 蛋白脑组织样蛋白 1 抗体

别名： FLJ25797; GTPase RhebL1; MGC34869; Ras homolog enriched in brain like 1; Ras homolog enriched in brain like 1 c; Ras homolog enriched in brain like protein 1; Ras homolog enriched in brain like-1 c; Ras homolog enriched in brain-like protein 1; REBL1_HUMAN; Rheb like protein 1; Rheb-like protein 1; Rheb2; Rhebl1; RhebL1c.

研究领域： 细胞生物 神经生物学 信号转导 表观遗传学 G 蛋白信号

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat, Cow, Rabbit,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 21kDa

细胞定位： 细胞浆 细胞膜

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human RHEBL1:101-183/183

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20 ° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed