

RAD51 相互作用蛋白 1 抗体

产品货号： mlR19725

英文名称： RAD51AP1

中文名称： RAD51 相互作用蛋白 1 抗体

别名： 2510006L10Rik; MGC156698; OTTMUSP00000028651; PIR51; R51A1_HUMAN; RAB22; RAD51 associated protein 1; RAD51 interacting protein; RAD51-associated protein 1; RAD51-interacting protein; Rad51ap1.

研究领域： 细胞生物 转录调节因子 表观遗传学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 38kDa

细胞定位： 细胞核

性状： Lyophilized or Liquid

浓度： 1mg/ml

免 疫 原： KLH conjugated synthetic peptide derived from human RAD51AP1:1-100/352

亚 型： IgG

纯化方法： affinity purified by Protein A

储 存 液： 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件： Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed： PubMed