

PRRC2C 蛋白抗体

产品货号： mlR19425

英文名称： PRRC2C

中文名称： PRRC2C 蛋白抗体

别名： 1810043M20Rik; 9630039I18Rik; A630006J20; BAT2 domain-containing protein 1; BAT2-iso; Bat2d; BAT2D1; BAT2L2; E130112L15Rik; HBV X-transactivated gene 2 protein; HBV XAg-transactivated protein 2; HLA-B-associated transcript 2-like 2; mKIAA1096; PRC2C_HUMAN; Proline-rich and coiled-coil-containing protein 2C; Protein PRRC2C; Prrc2c; Prrc3; RGD1566064; XTP2.

研究领域： 肿瘤 细胞生物 免疫学

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Dog, Pig, Cow, Horse, Rabbit, Zebrafish, Sheep, Rhesus monkey, Orangutan, Elephant

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分子量： 317kDa

细胞定位： 细胞浆

性状： Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human PRRC2C:1-100/2896

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed