

## PRA1 蛋白抗体

产品货号： mlR19980

英文名称： PRA1/PRAF1

中文名称： PRA1 蛋白抗体

别 名： PRA1 domain family 1; PRA1 family protein 1; PRAF1; PRAF1\_HUMAN; prenylated Rab acceptor 1; Prenylated Rab acceptor protein 1; Rab acceptor 1 (prenylated); RABAC1; YIP3.

研究领域： 细胞生物 信号转导 G 蛋白信号

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Mouse, Rat,

产品应用： ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500 （石蜡切片需做抗原修复）

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量： 21kDa

细胞定位： 细胞浆 细胞膜

性 状： Lyophilized or Liquid

浓 度： 1mg/ml

**免 疫 原：** KLH conjugated synthetic peptide derived from human PRA1:1-100/185 <Cytoplasmic>

**亚 型：** IgG

**纯化方法：** affinity purified by Protein A

**储 存 液：** 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**保存条件：** Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

**PubMed：** PubMed