

颌骨成釉细胞相关蛋白抗体

产品货号： mlR21131

英文名称： ODAM

中文名称： 颌骨成釉细胞相关蛋白抗体

别名： 2310011G06Rik; Apin; FLJ20513; ODAM; ODAM_HUMAN; Odontogenic ameloblast associated; Odontogenic ameloblast-associated protein.

研究领域： 细胞生物 发育生物学 信号转导 细胞周期蛋白

抗体来源： Rabbit

克隆类型： Polyclonal

交叉反应： Human, Chimpanzee,

产品应用： WB=1:500-2000 ELISA=1:500-1000 IHC-P=1:400-800 IHC-F=1:400-800 ICC=1:100-500 IF=1:100-500
(石蜡切片需做抗原修复)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

分 子 量 : 29kDa

细胞定位 : 细胞核 细胞浆 分泌型蛋白

性 状 : Lyophilized or Liquid

浓 度 : 1mg/ml

免 疫 原 : KLH conjugated synthetic peptide derived from human ODAM:181-279/279

亚 型 : IgG

纯化方法 : affinity purified by Protein A

储 存 液 : Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4

保存条件 : Store at -20 ° C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20° C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 ° C.

PubMed : PubMed

产品介绍 Function:

Tooth-associated epithelia protein that probably plays a role in odontogenesis, the complex process that results in the initiation and generation of the tooth. May be incorporated in the enamel matrix at the end of mineralization process.

Subcellular Location:

Secreted.

Post-translational modifications:

O-glycosylated.

Similarity:

Belongs to the ODAM family.

SWISS:

A1E959

Gene ID:

54959

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

产品图片

